otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

C € 0598

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Thin nitrile disposable glove with beaded cuff. The glove is powder free and ambidextrous.

Generally

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

Category

The glove is certified in category III in compliance with Medical Device Directive 93/42/EEC and Personal protection Equipment Regulation (EU) 2016/425. The glove has been tested in accordance with the standards EN 420:2003+A1:2009 (general requirements) and EN ISO 374-1:2016, EN 374-4:2013, EN ISO 374-5:2016 (chemicals and micro organisms). Notified body responsible for EU Type Examination and Module D conformity assessment: SGS FIMKO OY, P.O. Box 30 (Särkiniementie 3) 00211 Helsinki, Finland. Notified Body number: 0598.

EN ISO 374-1:2016/Type B

KPT

EN ISO 374-5:2016

VIRUS

EN 420: These gloves do not comply with EN 420 hand sizing. Fit for special purpose only.**EN ISO 374-1:2016/Type B****EN 374-4:2013**

	Chemical	Time*	Level	Degradation %
A	Methanol	>2	0	80.6
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A
I	Ethyl acetate	N/A	N/A	N/A
J	N-heptane	>13	1	33.2
K	Sodium hydroxide 40%	>480	6	8.9
L	Sulphuric acid 96%	>5	0	100.0
M	Nitric acid 65 %	N/A	N/A	N/A
N	Acetic acid 99 %	N/A	N/A	N/A
O	Ammonium hydroxide 25%	>10	1	80.8
P	Hydrogen peroxide 30 %	>60	3	N/A
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyd 37 %	>44	2	33.7

*)Breakthrough time (minutes)

**)Permeation rate: 7 µg/cm2/min

Performance Level	1	2	3	4	5	6
Minimum breakthrough time (mins)*	10	30	60	120	240	480

*) Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

EN374-4:2013

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

EN ISO 374-5:2016Resistance to Bacteria and Fungi = Pass
Resistance to Virus = PassThe glove is approved for food handling.
(Regulation (EC) No. 1935:2004)**Application / qualities**

Thin and flexible seamless glove with great dexterity. The glove is waterproof and suitable for handling of many different chemicals with short contact time. If the glove has been in contact with substances it must be discarded, when the breakthrough time is reached. Primarily used in food production, mounting, product control, laboratories, chemical industry and for general cleaning.

Treatment / storage

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging. When stored as recommended, properties of the glove will not be changed up to five years.

Warning

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimens.

Packaging100 gloves per dispenser.
10 dispensers in a carton made of recyclable cardboard.**SE • Bruksanvisning****Varumärke / art. nr.****Nilex PF**608070 Storlek 6½-7
608080 Storlek 7½-8
608090 Storlek 8½-9
608090 Storlek 9½-10
608110 Storlek 11**Beskrivning**

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otto schachner

Prinsessens Kvarter 2
7000 Fredericia, Denmark

EU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

C € 0598

otto schachner
member of ENEC Group

C	Acetonitril	N/A	N/A	N/A
D	dichlormethan	N/A	N/A	N/A
E	Carbondisulfid	N/A	N/A	N/A
F	Toluen	N/A	N/A	N/A
G	Diethylamin	N/A	N/A	N/A
H	tetrahydrofuran	N/A	N/A	N/A
I	Ethylacetat	N/A	N/A	N/A
J	N-heptan	>13	1	33.2
K	Natriumhydroxid, 40 %	>480	6	8.9
L	Svovlsyre, 96 %	>5	0	100.0
M	Salpetersyre 65 %	N/A	N/A	N/A
N	Eddikesyre 99 %	N/A	N/A	N/A
O	Ammoniumhydroxid 25%	>10	1	80.8
P	Brintoverilte 30 %	>60	3	N/A
S	Fluorsyre 40 %	N/A	N/A	N/A
T	Formaldehyd 37 %	>4	2	33.7

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