**GE • User instructions**

**Nitrile Disposable Gloves, 0.11-0.16 mm**

Blue, Art. No.	Black, Art. No.	Size
7800884050	7800884054	8/M
7800884051	7800884055	9/L
7800884052	7800884056	10/XL
7800884053	7800884057	11/XXL

**Description**

Nitrile single use examination glove with beaded cuff. The glove is powder free and ambidextrous.

**Quality:**

Weight: 5.5 g  
Thickness of palm: 0.11 - 0.12 mm  
Thickness of fingertips: 0.15 – 0.16 mm  
AQL 1.5

**General**

When donning the glove, please ensure that both glove and hand are clean, glove is free from defects that can hinder performance, glove size is right, and glove is properly fit on the contours and crotches of fingers. In case of contamination / perspiration, take of the glove, allow it to dry before wearing again and/or discard depending on the condition of glove.

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

**Category**

The glove is certified in category III in compliance with Medical Device Regulation 2017/745 and PPE Regulation (EU) 2016/425 on safety requirements for personal protective equipment. The gloves were tested in accordance with the regulation no. 10/2011/EU and EU regulation no. 1935/2004 on Contact with foodstuff – Please check the specifications at the end of this document.

The glove has been tested in accordance with the standards EN ISO 21420:2020 (general requirements), EN 455 "Medical gloves for single use": EN 455-1:2020 (freedom from holes), EN 455-2:2020 (physical properties), EN 455-3:2015 (biological evaluation), and EN 455-4:2009 (shelf-life determination). EN 16523-1:2015+A1:2018 (resistance to permeation of chemicals), and EN ISO 374-1:2016+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016 (chemicals and micro-organisms). Notified Body responsible for EU Type Examination and ongoing supervision (Module D): SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland (2777).

EN ISO 374-1:2016+A1:2018/Type B EN ISO 374-5:2016

KPT

VIRUS

**NITRIL-ENGANGSHANDSKER  
PUDDERFRI**

Farver: BLÅ og SORT

**EN ISO 374-1:2016+A1:2018/Type B**

	Chemical	Time*	Level	Degradation %
A	Methanol	N/A	N/A	N/A
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A
I	Ethyl acetate	N/A	N/A	N/A
J	N-heptane	N/A	N/A	N/A
K	Sodium hydroxide 40%	>480	6	-11.5
L	Sulphuric acid 96%	N/A	N/A	N/A
M	Nitric acid 65 %	N/A	N/A	N/A
N	Acetic acid 99 %	N/A	N/A	N/A
O	Ammonium hydroxide 25%	N/A	N/A	N/A
P	Hydrogen peroxide 30%	>30	2	-9.5
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyde 37 %	>60	3	7.4

\*Breakthrough time (minutes)

Permeation Performance level	1	2	3	4	5	6
Measured breakthrough time (mins)*	>10	>30	>60	>120	>240	>480

\*) Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

**EN ISO 374-4:2019**

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

**EN ISO 374-5:2016**

Resistance to Bacteria and Fungi = Pass  
Resistance to Virus = Pass

**Migration test**  
(EN 1186-1:2002)

**Application / qualities**

Thin and flexible seamless glove with great dexterity. The glove is waterproof and suitable for handling of many different chemicals with short contact time. If the glove has been in contact with substances it must be discarded, when the breakthrough time is reached. Primarily used in food production, mounting, product control, laboratories, chemical industry and for general cleaning.

**Treatment / storage**

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging.

**EN ISO 374-4:2019**

**Warning, EN ISO 374-1:2016+A1:2018:**

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.
- The glove does not contain any substances that are known to cause allergies.

**Packaging**

100 gloves per dispenser.  
10 dispensers in a carton made of recyclable cardboard.

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The gloves are in conformity with EU Regulation 1935/2004:

**Test Requested and Conclusion(s):**

No.	Test Sample	Standard and Requirement	Conclusion(s)
1)	Submitted sample	According to Regulation EU 1935/2004/EC & EU 10/2011& (EU) 2020/1245. - Sensory test - Overall Migration - Specific migration of Primary Aromatic Amine (PAA) - Specific Migration of Heavy metals - Specific migration of Nitrosamine and nitrosatable substances	PASS
2)	Submitted sample	German Food, Articles of Daily Use and Feed Code of September 1, 2005 (LFGB), Section 30 & 31 with amendments and BfR recommendation III and BfR Recommendation L&VI. & European Commission Regulation (EC) No 10/2011 with amendments. - Total butadiene	PASS

**Test Result(s):**

**1) According to Regulation EU 1935/2004/EC & EU 10/2011 & (EU) 2020/1245.**

**1.1) Sensory test**

Method: With reference to DIN 10955:2004.

No. of panelist: 6.

Test Item(s)	Test condition	Limit	Result	Conclusion
			1	
Sensorial examination odour (Point scale).	23°C, 24h	2.5	1	PASS
Sensorial examination taste (Point scale).	23°C, 24h Distilled water	2.5	0	PASS

Notes: Scale evaluation:

0= no aberration, neutral

1= very slight deterioration, barely perceptible

2= slight deterioration

3= significant deterioration

4= strong deterioration

**1.2) Overall migration**

Method: EN 1186-1:2002 & EN 1186-3:2002 & EN 1186-2:2002

Test condition: 40°C for 0.5h in 10% Ethanol / 3% Acetic acid / Olive oil.

Simulant Used(s)	Unit	Limit(s)	RL	Result(s)	Conclusion(s)
				1	
10% Ethanol	mg/dm <sup>2</sup>	10	3	N.D.	PASS
3% Acetic acid	mg/dm <sup>2</sup>	10	3	3.8	PASS
Olive oil	mg/dm <sup>2</sup>	10	3	N.D.	PASS

**Note:** 1. mg/dm<sup>2</sup> = milligram per square decimeter.

2. N.D. = Not Detected (< RL).

3. RL = Reporting Limit.

4. Solvents & test condition were selected based on Client's request.

The tested items **complied with** the requirements of German Food & Feed Acts of September 1, 2005 (LFGB), Section 30 and 31.

- Overall migration test

- Extractable formaldehyde

- Total lead and zinc content

- Color release

- Sensory test

As per client's request

- Total Cadmium content

**001 PASS**

**001 PASS**

**001 PASS**

**001 PASS**

**002 PASS**

**001 PASS**

**001 PASS**

**001 PASS**

**1. Overall Migration Test**

- Test method: In accordance with REGULATION (EU) No 10/2011 and its amendments on plastic materials and articles intended to come into contact with food.

- As specified in REGULATION (EU) No 10/2011 and its amendments; with reference to EN 1186: Part 4 (Test methods for overall migration into olive oil by cell)/ EN 1186: Part 5 (Test methods for overall migration into aqueous food stimulants by cell)

- Migration ratio (S/V): 10dm<sup>2</sup>/L

Simulant(s) Used	Test Condition	Result(s) [mg/dm <sup>2</sup> ]	Maximum Permissible Limit [mg/dm <sup>2</sup> ]
		001	
3% Acetic acid	40°C for 0.5 hour	4.80	10
10% Ethanol	40°C for 2 hours	4.20	10
Olive Oil	40°C for 2 hours	<0.500	10

**2. Extractable formaldehyde**

- Test method: For compliance with the Recommendation of the BfR "Kunststoffe im Lebensmittelverkehr" Part XXI. Commodities based on Natural and Synthetic Rubber

- With reference to Section 2.7.1 of methods for the "Testing of commodities made of rubber"

- Test condition: 3% Acetic acid, 40°C for 0.5 hour

- Migration ratio (S/V): 8dm<sup>2</sup>/L

Test Item(s)	Result(s) [µg/ml]	Maximum Permissible Limit [µg/ml]
	001	
Extractable Formaldehyde	<0.10	3

**3. Total lead and zinc content**

- Test method: Acid digestion, then followed by ICP-OES

Test Item(s)	Result(s) [%]	Maximum Permissible Limit [%]
	001	
Lead content	<0.001	0.003
Zinc content	0.525	3.0

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**5. Sensory test**

- Test method: With reference to DIN 10955.

- The submitted sample was immersed in distilled water at 40°C for 2 hours. After this treatment treated water was examined by panels with regard to any divergence in smell and taste.

Sample(s)	Testing Parameter	Grading result(s)	Recommended level
002	Transfer of taste	0.5	<3
	Transfer of smell	0.5	<3

Note:

1. Available grading are listed as follow:

Grading 0: No perceptible taste/smell deviation

1: Just perceptible taste/smell deviation

2: Weak taste/smell deviation

3: Clear taste/smell deviation

4: Strong taste/smell deviation

**6. Total Cadmium content**

Test Item(s)	Result(s) [mg/kg]	Maximum Permissible Limit* [mg/kg]
	001	
Total Cadmium Content	<2	2

Testdata re. Contact with Foodstuff: Regulation No. 10/2011/EU, EN1186, etc.: