CleanSpace[®]

R E S P I R A T O R S

USER INSTRUCTIONS

www.cleanspacetechnology.com

■ CleanSpace HALO CS3000

CleanSpace Technology Pty Ltd
ABN 24 146 453 554

T: +61 2 8436 4000

GF, 16-18 Carlotta St
Artarmon NSW 2064 Australia

sales@cleanspacetech.com
nology.com

CleanSpace®

R E S P I R A T O R S

CONTENT - LANGUAGES

English.....	2-30
Croatian.....	31-60
Czech.....	61-90
Danish.....	91-120
Dutch.....	121-153
Estonian.....	154-183
Finnish.....	184-213
French.....	214-243
German.....	244-276
Hungarian.....	277-308
Italian.....	309-341
Latvian.....	342-372
Lithuanian.....	373-403
Norwegian.....	404-433
Polish.....	434-464
Portuguese.....	465-495
Romanian.....	496-526
Russian.....	527-559
Serbian.....	560-591
Slovak.....	592-622
Bulgarian.....	623-652
Spanish.....	653-685
Swedish.....	686-715
Turkish.....	716-745
Ukrainian.....	746-778

1. Contents

1.	Contents	2
2.	Foreword	2
3.	System Description	3
4.	List of Warnings within these User Instructions 	3
5.	Special or Critical User Instructions	4
6.	Limitations of Use	4
7.	System Components	6
8.	Controls and Indicators	7
9.	Mask Fitting	8
10.	Using CleanSpace HALO With Half Mask	9
11.	Using CleanSpace HALO With Full Face Mask	17
12.	Entering and Exiting the Contaminated Area	22
13.	Working in Noisy Environments	23
14.	Fitting & Changing the Filter	23
15.	Battery Information	24
16.	Cleaning	26
17.	Periodic Maintenance & Checking	28
18.	Approvals	28
19.	Specifications	28
20.	Product and Accessory information	29
21.	Product Warranty	29

Important: Before use, the wearer must read and understand these user instructions. Keep these user instructions for reference.

Read these instructions in conjunction with the appropriate CleanSpace™ filters and accessories instructions.

WARNING
*This product is part of a system that helps protect against certain airborne contaminants. **Misuse may result in sickness or death.** For proper use, see your Occupational Health Specialist, these user instructions or contact CleanSpace Customer Service on 1 888 804 0038*

2. Foreword

Read all instructions and warnings before using this device. Keep these user instructions for reference. If you have questions regarding this system contact CleanSpace Customer Service sales@cleanspacetechnology.com or visit our website: www.cleanspacehealth.com

3. System Description

The components of the CleanSpace HALO system are shown in Section 7. The respirator draws ambient air through the high efficiency filter and supplies filtered air to the wearer via the mask. CleanSpace HALO is a breath responsive respirator which means that it blows harder during inhalation and more softly during exhalation. The device continuously adjusts fan speed to maintain positive pressure in the mask. Should the filter become blocked, an audible alarm will sound and the “Filter” light on the keypad will be illuminated. If the battery voltage falls below the minimum required to sustain the designed flow rate, an audible alarm will sound.

WARNING

*Properly selected, used and maintained respirators help to protect against certain airborne contaminants. It is essential to follow all instructions and government regulations on the use of this product, including wearing the complete respirator system during all times of exposure, in order for the product to help protect the wearer. **Misuse of respirators may result in overexposure to contaminants and lead to sickness or death.** For proper use, see your Occupational Health Specialist, these user instructions or contact CleanSpace Customer Service sales@cleanspacetechnology.com or visit our website: www.cleanspacehealth.com*

4. List of Warnings within these User Instructions

WARNING

*This product is part of a system that helps protect against certain airborne contaminants. **Misuse may result in sickness or death.** For proper use, see your Occupational Health Specialist, these user instructions or contact CleanSpace Customer Service on sales@cleanspacetechnology.com or visit our website: www.cleanspacehealth.com*

WARNING

*Properly selected, used and maintained respirators help to protect against certain airborne contaminants by reducing concentrations in the wearer's breathing zone. It is essential to follow all instructions and government regulations on the use of this product, including wearing the complete respirator system during all times of exposure, in order for the product to help protect the wearer. **Misuse of respirators may result in overexposure to contaminants and lead to sickness or death.** For proper use, see your Occupational Health Specialist, these user instructions or contact CleanSpace Customer Service on sales@cleanspacetechnology.com or visit our website: www.cleanspacehealth.com*

**CleanSpace HALO is a tight fitting respirator and National standards recommend the user be clean shaven.
IF A SATISFACTORY FIT FACTOR CANNOT BE ACHIEVED WITH ANY OF THE MASKS, CLEANSPEACE HALO MUST NOT BE USED**

**DO NOT USE COMPRESSED AIR OR A BRUSH TO CLEAN THE FILTER!
FILTERS ARE DAMAGED BY THE USE OF COMPRESSED AIR OR BY BRUSHING.**

**Misuse of CleanSpace Respirators by cleaning the filter may result in overexposure to contaminants and lead to sickness or death.
IF A SATISFACTORY FIT CANNOT BE ACHIEVED, DO NOT ENTER THE CONTAMINATED ZONE**

Be sure to remove the Seal Check Cap before entering the contaminated area. The Cap blocks the exhalation valve, making it more difficult for your exhaled air to be expelled from the mask. Failure to remove the cap will lead to build-up of Carbon Dioxide in the mask, and may result in headache or dizziness. Never leave the Cap in place for more than 2 minutes!

DO NOT REMOVE THE RESPIRATOR UNTIL YOU HAVE VACATED THE CONTAMINATED AREA unless you have immediate and sudden onset health reasons to do so (for instance you are experiencing dizziness and believe removing the respirator while you leave the contaminated area may help)

IT IS ESSENTIAL THAT THE CORRECT FILTER TYPE IS SELECTED FOR THE CHOSEN APPLICATION.

If the blocked filter alarm is triggered (2 beeps, repeated every second, red LED flashes), leave the contaminated area immediately and replace the filter. Operating the respirator after the blocked filter alarm has sounded can cause the flow to fall below the manufacturer's minimum designed flow, which may result in overexposure to contaminants and lead to sickness or death.

**DO NOT USE WATER, COMPRESSED AIR OR ANY OTHER MECHANICAL METHOD TO CLEAN THE FILTER!
FILTERS ARE EASILY DAMAGED BY THE USE OF WATER, COMPRESSED AIR OR BY BRUSHING.**

Misuse of CleanSpace Respirators by cleaning the filter may result in overexposure to contaminants and lead to sickness or death.

Use only NIOSH or CE Mark-approved CleanSpace™ filters with the CleanSpace HALO™ respirator. Certification or approvals are void if used with other filters other than CleanSpace. Use of other filters may result in overexposure to contaminants and lead to sickness or death.

CleanSpace HALO is not intrinsically safe. Do not use in flammable or explosive atmospheres. Doing so may result in injury or death.

Always correctly use and maintain the internal lithium ion battery packs. Failure to do so may result in fire or explosion or could adversely affect respirator performance and result in injury, sickness or death.

- Do not charge the internal battery with unapproved chargers or in enclosed cabinets without ventilation or near flammable liquids or gasses, or near sources of high heat
- Do not immerse the device in water
- Do not use, charge or store the device outside the recommended temperature limits

If the battery alarm sounds (3 beeps, repeated every second), leave the contaminated area immediately and re-charge the battery. Operating the respirator after the low battery alarm has sounded can cause the flow to fall below the manufacturer's minimum designed flow, which may result in overexposure to contaminants and lead to sickness or death.

The Flow Test is not a battery charge test. Three LEDs means that, as of the moment it is tested, the unit can deliver high flow. It does not mean the battery is full. You must check the battery condition separately. See Section 16.

IN THE EXTREMELY RARE CIRCUMSTANCE THAT THE BATTERY IS DAMAGED AND ELECTROLYTE COMES IN CONTACT WITH EYES, FLUSH WITH WATER IMMEDIATELY AND SEEK URGENT MEDICAL ATTENTION

WARNING: NEVER DRY THE MASK OR EXHALATION VALVE WITH A CLOTH THAT MAY LEAVE BEHIND LINT.

Lint contamination of the exhalation valve may cause it to leak, resulting in overexposure to contaminants and lead to sickness or death.

WARNING: The Respirator contains a battery, sensitive electronics and a motor. NEVER immerse it in water or use anything wetter than a damp cloth to clean it.

WARNING: NEVER CLEAN THE FILTER WITH WATER, COMPRESSED AIR OR A BRUSH. FILTERS ARE DAMAGED BY THE USE OF WATER COMPRESSED AIR OR BY BRUSHING.

Misuse of CleanSpace HALO by cleaning the filter may result in overexposure to contaminants and lead to sickness or death.

You must recalibrate the internal pressure sensor any time that your CleanSpace Respirator is exposed to changes in temperature of more than 20°C

5. Special or Critical User Instructions

This respirator contains a system for synchronising with your breathing and regulating mask pressure. **This system requires re-calibration whenever it experiences a change in temperature of more than 20°C. It is best practice to also re-calibrate if the unit has been in storage, particularly if the storage temperature is not known.** For the most accurate battery charge indication, the respirator must be plugged into its charger. To be sure it is fully charged, plug in the charger and wait for **all 3 LEDs to light solidly without flashing**. If the 3rd LED is flashing rapidly, the battery is 95% charged.

6. Limitations of Use

Use this Respirator strictly in accordance with all instructions in these user instructions. Never modify or alter this product.

- **Do not remove the respirator until you have left the contaminated area**, unless you have immediately and sudden onset health reasons to do so (for example you are experiencing dizziness and believe removing the respirator while you leave the contaminated area may help).
- CleanSpace HALO is a tight fitting respirator and National standards recommend the user be clean shaven. If a satisfactory fit factor cannot be achieved with any of the masks, CleanSpace Respirators must not be used
- This Respirator WITH THE HALF MASK DOES NOT PROTECT YOUR EYES. In conditions that may damage or irritate eyes, use protective eyewear.
- Only use your respirator with the parts and accessories listed on the Approval Label.
- Do not use the respirator unless it is powered and running normally.
- Do not use the respirator while it is being charged via the Battery Charger
- Do not use in airborne contaminant concentrations above those specified in your national regulations
- Do not use for respiratory protection against unknown atmospheric contaminants or when concentrations of contaminants are unknown or immediately dangerous to life or health (IDLH).
- Do not use in oxygen deficient or oxygen enriched atmospheres. Do not use in flammable or explosive environments.
- Only for use by trained personnel.

- Filters need to be changed regularly. The frequency of change depends on use and the concentration of contaminants in the atmosphere.
- Do not use for escape purposes. National regulations may impose specific limitations on the use of filters depending on the filter class and the facemask used.

LEAVE THE CONTAMINATED AREA IMMEDIATELY IF:

- Respirator warning lights and/or sounds activate for low battery or blocked filter
- Any part of the respirator is damaged
- Air flow into the mask decreases or stops
- Breathing becomes difficult or increased resistance occurs
- You feel dizzy or your airway is irritated
- You can taste or smell contaminants

CleanSpace respirators are suitable for use in the following atmospheric conditions:

- Operating Temperature Range: -10°C to 45°C
- Relative humidity: 0 to 90% non-condensing

The respirator will stop functioning if its internal temperature rises above 60°C or falls below -10°C

If the respirator has been used in an area that has caused it to become contaminated with a substance requiring special decontamination procedures it should be placed in a suitable container and sealed until it can be decontaminated.

Do not disassemble the respirator case. There are no user serviceable parts inside.

Failure to follow all instructions on the use of this product, and/or failure to use the respirator during times of exposure, may lead to adverse effects on the wearer's health and may render the warranty void.

7. System Components


Figure 1 - System Components

PARTS OF THE RESPIRATOR

Throughout this manual reference is made to various commonly-used components and features of the Respirator. Familiarise yourself with these parts before reading the rest of the manual. See Figure 2.


Figure 2 - Parts of the Respirator

8. Controls and Indicators

POWER BUTTON

This button is used to switch between the three (3) operating Modes: “On”, “Standby” and “Off” Mode.

i) **Standby Mode** is when the green battery indicator lights are on, the motor is not running and there is no airflow to the mask. The respirator will automatically switch to Standby Mode within ten (10) seconds of the user taking it off. If the respirator is in Off Mode it will switch to Standby Mode when the Power button is pressed.

ii) **On Mode** is when the respirator is being worn, the motor is running and there is air flow to the mask. The motor will start (called On Mode) when the respirator detects a change in pressure in the mask triggered by your breathing. In On Mode you should hear the motor running and feel the airflow on your face. You can also switch to On Mode (start the motor) from Standby Mode by pressing the Power button once.

iii) **Off Mode** is when the green battery indicator lights are off and the motor is not running and there is no air flow to the mask. The respirator automatically switches into Off Mode three (3) minutes after it has been removed from the user's face and the sensors detect that there is no breathing.

To conserve battery life, the respirator is designed to automatically switch from On Mode to Standby then to Off Mode when not being worn.

Important: In Standby Mode, the CleanSpace Respirator switches to On Mode as the wearer starts to breathe. The respirator can only switch into On Mode from Standby Mode. If the respirator is in OFF mode, the motor will not run.

FLOW TEST BUTTON

This button is used to check that the respirator is able to deliver its minimum designed flow. Pressing it once, when the CleanSpace is in Standby Mode (not ON Mode), starts the flow test which lasts about 5 seconds. **See Section 10 for instructions on running the flow test.**

BATTERY INDICATOR LIGHTS

Your respirator is equipped with an indicator of battery charge. There are three (3) battery indicator lights. With the charger plugged in and three green lights lit, the battery is fully charged and typically has up to eight (8) hours of operating time. The respirator should be fully charged before use.

If you need to be sure the battery is 100% charged plug in the charger. Even if the battery is fully charged, the 3rd LED will flash for at least three minutes while the respirator checks its condition. Once **all 3 LEDs light solidly without flashing the battery is at 100%**

For how to assess the level of charge, see Section 14.

When the battery approaches a level at which it would not be able to supply the Manufacturer's Minimum Design Flow, an alarm sounds (3 beeps, repeated once per second). All green battery lights are extinguished. **If the low battery alarm sounds you must leave the contaminated area immediately and recharge the battery.**

Operating time is affected by mask seal, work rate and other factors. The operating times quoted above are average durations at moderate work rates at sea level. Actual operating times may vary widely from average durations.

FILTER WARNING ALARM

CleanSpace Respirators have a Filter Warning Alarm, which is triggered when the filter is blocked. **If the Filter Warning Alarm sounds** (two beeps, repeated once per second) or the Filter Warning Alarm light comes on, **you must move out of the contaminated area, and change the filter.**

9. Mask Fitting

Before using your respirator you must determine the right mask for your face and know how to adjust the system to achieve a good fit.

It is vital that your mask is the right size for your face and fits properly. **Mask fitting must be carried out by a specialist / designated mask fitter.**


To carry out a quantitative fit test your specialist mask fitter will require **PAF-0025 PortaCount Adaptor** and the accompanying **S005-7174 Quantitative Fit Test Instructions**. The PortaCount Adaptor is an accessory that fits between the mask and the exhalation valve and allows the air in the mask to be sampled. It comes with detailed instructions on how to carry out a quantitative fit test.

CleanSpace HALO is a tight fitting respirator and National standards recommend the user be clean shaven. IF A SATISFACTORY FIT FACTOR CANNOT BE ACHIEVED WITH ANY OF THE MASKS, CLEANSPEACE HALO MUST NOT BE USED

MASK SELECTION GUIDELINES

No set of guidelines can ensure that you have the right size mask for your face. You must confirm the fit with a quantitative fit test. But the table below will help your mask fitting specialist to select the mask most likely to fit you, and may therefore save time.

Once you have achieved an acceptable mask fit and confirmed it with a quantitative fit test, record your information in the table below.

Name	
Mask Size	
Neck Support	Fitted / Not Fitted (Circle one)
Date of Quantitative Fit Test	
Fit Factor on Test	

IF A SATISFACTORY FIT FACTOR CANNOT BE ACHIEVED WITH ANY OF THE MASKS, CLEANSPACE RESPIRATORS MUST NOT BE USED

10. Using CleanSpace HALO With Half Mask

Complete the following five steps each time you use your respirator.

STEP 1 - INSPECT

Before each entry into a contaminated area, the following inspections must be performed:

- Visually check the all parts including the Respirator, mask, harness and filter. If parts are missing or damaged replace only with approved parts before proceeding.
- Check the RESPIRATOR** for cracks, holes or other damage or missing parts. Do not use the device if there is any damage or misuse. Check both bellows for splits or holes. Check that the bellows have not become distorted so as to partially or fully close the air path to the mask. Check and fit the neck support (if required). Check the battery is fully charged by pressing the Power button. All three (3) battery Green LEDs must light up. See Section 14.

- **Check the FILTER** carefully. The foam seal must be clean and free from damage of any kind. Examine the visible internal surfaces for any sign that dust has penetrated the media. If dust is found, do not use the filter. The body of the filter must be free of cracks or signs of damage. If any sign of impact or scratching is found, discard the filter. Fit the filter to the Respirator (see Changing the Filter, below)

**DO NOT USE COMPRESSED AIR OR A BRUSH TO CLEAN THE FILTER!
FILTERS ARE DAMAGED BY THE USE OF COMPRESSED AIR OR BY BRUSHING.
Misuse of CleanSpace Respirator by cleaning the filter may result in overexposure to
contaminants and lead to sickness or death.**

- **Check the MASK** to ensure that there are no cracks, tears or dirt; check the mask is not distorted. Check the mask exhalation valve for damage or dirt build up. If it is dirty, remove the exhalation valve cover. Remove any dirt, hairs or anything that could affect the seal of the valve against its seat. Check that the valve seat is clean. Reinstall the valve cover. If the valve is damaged, replace the mask.
- **Check the HEAD HARNESS** is intact and has good elasticity and can be fitted to the attachments on the mask. It must be adjusted to support the mask to seal comfortably to the face. DO NOT overtighten.

STEP 2 – CALIBRATE

This respirator contains a system for synchronising with your breathing and regulating mask pressure. **The system should be re-calibrated if it experiences a significant change in temperature of more than 20°C. It is best practice to re-calibrate if the unit has been in storage when the storage temperature is unknown.**

Steps to re-calibrate:

1. Remove the filter from the respirator.
2. With CleanSpace Respirator in Standby Mode (ie one or more green LEDs lit), press the “Power” and “Test Flow” buttons on the keypad at the same time.
3. The red and blue LEDs will both light, and the green battery LEDs will cycle.
4. Hold the respirator still until the red and blue LEDs turn off (10 – 15 seconds). After this, calibration is complete.
5. At the end of the calibration process the green LEDs return to indicating battery charge status.

STEP 3 – TEST FLOW RATE

This test checks that the Respirator is able to deliver the Manufacturer's Minimum Design Flow of 120 litres/minute.

1. Remove the mask from the respirator. Place the respirator flat on a table or other support.
2. Press and release the button marked “Flow Test”
3. The respirator will automatically run the Flow Test. Note: During the test, the motor will run fast and air will blow from the left hand bellows.
4. After 2 seconds the respirator reports the Test result using the LEDs on the keypad. Use the table below to interpret the LEDs.

LIGHTS	Meaning
3 LEDs:	PASS (Excellent: flow >180 l/min)**
2 LEDs:	PASS (Good)**
1 LED:	PASS (Acceptable)
ALL LEDs FLASH	FAIL (Flow <120 l/min) Do not use the respirator until a new filter has been fitted and / or battery fully charged and the test has been repeated with a PASS result. Reset the respirator by pressing the Power button. Fully charge the battery and / or replace the filter. Repeat the flow test. If the respirator fails the test filter is new and battery fully charged but, contact CleanSpace Technology and do not use until it has been evaluated

****WARNING:** The Test Flow Rate is not a battery charge test. Three LEDs means that, as of the moment it is tested, the unit can deliver high flow. It does not mean the battery is full. You must check the battery charge separately. See Section 16.

STEP 4 - DON THE RESPIRATOR

NOTE – please read STEP 4 all the way through before beginning to don the respirator

Locate your respirator, mask, neck support and harness. **Make sure the mask is the size that you used for your last successful quantitative fit test.** This size should be recorded in this manual (see Section 9). **ENSURE A CLEAN FILTER IS FITTED TO THE POWER UNIT.**

- **Fit a CleanSpace mask to the respirator**

With the peak of the mask upright locate the right hand mask AirClip. Locate the right hand respirator AirClip. It is on the same side of the respirator as the Keypad.

Figure 3 - Right hand mask AirClip




Figure 4 - Right hand *respirator* AirClip

**ENSURE A CLEAN FILTER IS
FITTED TO THE POWER UNIT.
(See Section 13)**

*Join the AirClips, so that the mask is attached to the respirator.
See Figure 5.*

Figure 5 - Mask on Blower (Note - connection is on right hand side)

Check that the mask is the right way up! The pointed nose section of the mask must face up, in the same direction as the keypad. See Figure 6.


Figure 6 - Nose section of mask and keypad must both face up

*Leave the other mask arm and bellow free until you are fitting the Respirator.
Familiarise yourself with the Adjust Buttons and the Mask Release Buttons. See Figure 7*


Figure 7 - Mask Release Button and Adjust Button

- **Fit the Harness**

Holding the harness find the two keyhole clips at the end of each elastic strap. And the marker, indicating the front of the harness. Fit the harness clips to the buttons on the inside of each Mask Airclip. Ensure that the harness strap with the FRONT marker is towards the front of the mask and that there are no twists in the straps. See Figure 8. And Figure 9.


Figure 8 - Harness Fitting The two straps of the harness fit to the inside of the AirClips on the mask.


Figure 9 – Adjust Harness length By moving cleats up or down strap. **ENSURE THE MARKER IS LOCATED AT THE FRONT OF THE HEAD HARNESS.**

- **Fit a Neck Support**

Align the neck support with the respirator, the two arms with the rounded heads should be at the top. Locate the bottom arms of the neck support in the grooves at the bottom of the unit and slide neck support upwards until you hear a 'click'.



Figure 10 - Fitting neck support

- **Loosen both bellows to their widest opening**

Press the adjust button on the respirator and loosen the bellows on both sides so they are fully extended. See Figure 11.

Loosen/Adjust

Figure 11 - Adjust Button

- **Place the Respirator in Standby Mode**

Press the Power button once. The Green Battery LED lights will light up on the respirator. The motor will not come on and the air will not flow in Standby Mode.

Power button on keypad

- **Don the Respirator and switch to On Mode**

Place the respirator behind your neck with the bellows and mask resting down one side,. The harness should be hanging in front of mask (Figure 12 – 14). Fasten the AirClips on your left-hand side so both sides of the mask are connected to the Respirator. Pull the mask up to sit comfortably on your face (Figure 15). **Start to breathe to trigger the system into On Mode for the air to flow.** Pull the harness back and onto your head (Figure 16). Adjust harness and settle the mask and power unit so both are comfortable. Adjust the bellows and harness to achieve a comfortable fit and seal on the face. **DO NOT OVER TIGHTEN THE SYSTEM** (Figure 17)

Figure 12 – Blower behind neck


Figure 13 – Fasten AirClip

Figure 14 – Device fastened around neck

Figure 15 – Pull mask up and onto face

Figure 16 – Pull harness back and onto head

Figure 17 - Settle harness and adjust


ON MODE: From Standby Mode, when the mask is placed on the face, the wearers breathing will trigger the system to start (On Mode) and the air to flow. Breathe normally. The wearer should feel fresh air on the face and the motor adjusting to the breathing patterns.

If the motor does not start and air does not flow, check the system was in Standby Mode (with the Battery green LEDS on). If the system is Off, (ie not in Standby Mode and no green LEDS), remove the respirator, press the Power button to activate Standby Mode and don again as described Figures 12 – 17 above.

- **Check the Harness**

The Harness should take a little of the weight of the mask and support the respirator so that the respirator system sits level on your head. If necessary adjust the harness. This is done by shortening or lengthening the front straps of the harness, by moving the cleats up or down the straps. Remember to shorten or lengthen each side by the same amount.

Figure 18 - Correct Position for CleanSpace HALO


- **Adjusting the mask on your face**

Adjust the mask for comfort and until you feel no leaks between the mask cushion and your face.

To tighten the system, place one hand against the back of the respirator and the other over the front of the mask. Squeeze the respirator forwards and the mask back. You will hear a series of clicks as the mask is tightened.

To Loosen the system, press the adjust buttons on either side of the bellows, and pull the mask and respirator away from each other.


Figure 19 – Adjusting the Mask

STEP 5 – CHECK MASK SEAL

It is recommended to do a fit check every time you wear a CleanSpace Respirator. To conduct a check on the mask seal, use PAF-1009, the red Seal Check Cap, and fit the Cap over the exhalation valve in the mask. The Cap will snap into place. See Figure 20. After the cap is fitted and you have checked no air is flowing from the exhalation valve, don the mask and breathe normally. Using your fingers, feel around the edge of the mask for leaks which can be felt as a cool air flow over your finger. Wetting your fingers makes it easier to feel leaks. For greater sensitivity breathe out firmly to raise the mask pressure. It is recommended that during the Seal Check you move your head by tilting the head up and down and side to side. Adjust the system if leaks are detected. Once you can feel no leaks from the mask, the seal check is complete. While conducting the Mask Seal Check, never leave the cap in place for more than 2 minutes while the respirator is donned.

Figure 20 - Fitting Seal Check Cap

IF A SATISFACTORY FIT CANNOT BE ACHIEVED, DO NOT ENTER THE CONTAMINATED ZONE

Before entering the contaminated areas, remove the Seal Check Cap by gently **twisting** the handle so that one side of the cap comes loose from the exhalation valve. Be careful not to dislodge the exhalation valve cover. See Figure 21.


Figure 21 - Removing Seal Check Cap

Be sure to remove the Seal Check Cap before entering the contaminated area. The Cap blocks the exhalation valve, making it difficult for exhaled air to be expelled from the mask. Failure to remove the cap will lead to build-up of Carbon Dioxide in the mask, and may result in headache or dizziness. Never leave the Cap in place for more than 2 minutes!

11. Using CleanSpace HALO With Full Face Mask

Complete the following five (5) steps each time you use your respirator with the Full Face Mask.

STEP 1 - INSPECT

Before entry into a contaminated area, the following inspections must be performed.

- Visually check the all parts including the Respirator, mask, harness and filter. If parts are missing or damaged replace only with approved parts before proceeding.
- **Check the RESPIRATOR** for cracks, holes or other damage or missing parts. Do not use the device if there is any damage or misuse. Check both bellows for splits or holes. Check that the bellows have not become distorted so as to partially or fully close the air path to the mask. Check and fit the neck support (if required). Check the battery is fully charged by pressing the Power button. All three (3) battery Green LEDs must light up. See Section 14.
- **Check the FILTER** carefully. The foam seal must be clean and free from damage of any kind. Examine the visible internal surfaces for any sign that dust has penetrated the media. If dust is found, do not use the filter. The body of the filter must be free of cracks or signs of damage. If any sign of impact or scratching is found, discard the filter. Fit the filter to the Respirator (see Changing the Filter, below)

**DO NOT USE COMPRESSED AIR OR A BRUSH TO CLEAN THE FILTER!
FILTERS ARE DAMAGED BY THE USE OF COMPRESSED AIR OR BY BRUSHING.
Misuse of CleanSpace Respirator by cleaning the filter may result in overexposure to
contaminants and lead to sickness or death.**

- **Check the FULL FCE MASK** to ensure that there are no cracks, tears or dirt; Check the mask is not distorted. Check the mask exhalation valve for damage or dirt build up. If it is dirty, remove the exhalation valve cover. Remove any dirt, hairs or anything that could affect the seal of the valve against its seat. Check that the valve seat is clean. Reinstall the valve cover. If the valve is damaged, replace the mask.
- **Check the HEAD HARNESS** is intact and has good elasticity and the fastening clips can be tightened and loosened. The harness must be adjusted to support the mask to seal comfortably to the face. **DO NOT** overtighten. Finally check the airclips can be clicked/unclicked into the Respirator.

STEP 2 – CALIBRATE

This respirator contains a system for synchronising with your breathing and regulating mask pressure. **The system should be re-calibrated if it experiences a significant change in temperature of more than 20°C. It is best practice to re-calibrate if the unit has been in storage when the storage temperature is unknown.**

Steps to re-calibrate:

- Remove the filter from the respirator.
- With CleanSpace Respirator in Standby Mode (ie one or more green LEDs lit), press the “Power” and “Test Flow” buttons on the keypad at the same time.
- The red and blue LEDs will both light, and the green battery LEDs will cycle.
- Hold the respirator still until the red and blue LEDs turn off (10 – 15 seconds). After this, calibration is complete.
- At the end of the calibration process the green LEDs return to indicating battery charge status.

STEP 3 – TEST FLOW RATE

This test checks that the Respirator is able to deliver the Manufacturer’s Minimum Design Flow of 120 litres/minute.

- Remove the mask from the respirator. Place the respirator flat on a table or other support.
- Press and release the button marked “Flow Test”
- The respirator will automatically run the Flow Test. Note: During the test, the motor will run fast and air will blow from the left hand bellows.
- After 2 seconds the respirator reports the Test result using the LEDs on the keypad. Use the table below to interpret the LEDs.

LIGHTS	Meaning
3 LEDs:	PASS (Excellent: flow >180 l/min)**
2 LEDs:	PASS (Good)**
1 LED:	PASS (Acceptable)
ALL LEDs FLASH	FAIL (Flow <120 l/min) Do not use the respirator until a new filter has been fitted and / or battery fully charged and the test has been repeated with a PASS result. Reset the respirator by pressing the Power button. Fully charge the battery and / or replace the filter. Repeat the flow test. If the respirator fails the test filter is new and battery fully charged but, contact CleanSpace Technology and do not use until it has been evaluated

**** WARNING:** The Test Flow Rate is not a battery charge test. Three LEDs means that, as of the moment it is tested, the unit can deliver high flow. It does not mean the battery is full. You must check the battery charge separately. See Section 16.

STEP 4 - DON THE RESPIRATOR WITH FULL FACE MASK

Donning is best done with a partner that can help check adjustments and fit!

If you have long hair, it is recommended you tie the hair back so that it does not interfere with the Full face harness and the mask seal on the face. To don the Full Face Mask, follow the steps below.

Step 1: Loosen all five straps on the mask harness to the fullest extent. Hold the mask in one hand while you use the other hand to pull the harness back and away from the mask.

Step 2: Place your chin in the cup of the mask face seal And Pull the harness over your head	
Step 3: Adjust the top strap so that the hanger for the rear strap sits an inch above the ears	
Step 3: Settle your face into the mask face seal. Gently tighten each harness strap in turn, starting with the bottom straps. As you pull each strap, use your other hand to steady the mask on your face. DO NOT OVER TIGHTEN. Continue to adjust the straps until the mask face seal presses evenly on the face. <u>Check that the mask seal does not cross the hairline</u> Ensure the mask is adjusted so the inner mask sits comfortably around your nose and mouth.	

CHECK THAT THE MASK SEAL DOES NOT CROSS YOUR HAIRLINE. Check all the way around the mask seal, paying particular attention to your forehead and temples. **The mask seal must not cross the hairline.**

IF YOU CANNOT ADJUST THE MASK TO AVOID YOUR HAIRLINE, IT IS RECOMMENDED THE MASK SHOULD NOT BE WORN

STEP 5 – MASK SEAL CHECK

It is recommend to do a negative pressure seal check on the mask.

Mask Seal Check:
Using your thumbs, **block the airclip** air inlets on both sides of the mask. **Breathe in sharply.** If there is a good mask seal, you should have difficulty drawing air as you inhale. Listen for squeaking or whistling noises which indicate air is leaking past the seal. **Hold your breath for 10 seconds.** If the mask seal slowly recovers (mask moves away from your face) there is a leak. **Readjust the mask fit and repeat the negative pressure seal check.** **Unblock the air inlets and resume normal breathing.**



IF A SATISFACTORY FIT CANNOT BE ACHIEVED, DO NOT ENTER THE CONTAMINATED ZONE

STEP 6 – DON THE RESPIRATOR

NOTE: Read STEP 4 all the way through before beginning to don the CleanSpace Respirator



FITTING THE NECK SUPPORT:




With the Full Face Mask donned, locate your respirator and neck support. There are two sizes of Neck Support (Small and Large). **Make sure neck support is of the size that you used for your last successful quantitative fit test.** Your size should be recorded in this manual. Fit the neck support to the Respirator. The neck support has keyhole openings that snap over buttons on the respirator case. Place the neck support in position against the buttons and then press firmly towards the respirator to engage the buttons to snap into the keyholes.

FAMILIARISE YOURSELF WITH THE RESPIRATOR

1. Find the adjust Buttons and the Mask Release Buttons.
2. Loosen both bellows to their widest opening

Place the Respirator in Standby Mode Press the Power button once and the Green LED Battery Lights will come on.	
Don the respirator. Place the respirator behind your head, resting on your shoulders. Note: The keypad on the respirator must face up. Grasp the left-hand mask AirClip in one hand and the left-hand respirator AirClip in the other. Connect and click into place. Breathe to trigger the air flow.	
With breathing, the CleanSpace Respirator will start. When air is flowing into the mask, join the AirClips on the right-hand side to close the system. IMPORTANT: If the motor does not start, check the system is in Standby Mode (green lights on). Take the respirator off, press Power button to activate Standby Mode and don again as described.	
NOTE! If you cannot connect the air clips on both sides because the respirator pulls up tight against the back of your neck, swap to a neck support suited to a larger head size. See above. Then start the respirator donning process again.	
With the Full Face Mask, The respirator can rest loosely or be carried close to the neck. By adjusting the bellows on either side of the respirator you can change the respirators position to suit your requirements. Jogging movements: it may be more comfortable for the unit to sit snugly against the neck. Looking up: Wear the respirator away from your neck.	
To adjust the respirator, use the bellow/adjust system on the sides. To tighten: Steady the mask and squeeze the respirator toward the mask toward the back of the neck. You will hear clicking as the adjust mechanism moves to the new position. To loosen: Press the adjust buttons on each side and pull backwards on the respirator to move it away from the back of the neck.	Adjust both sides evenly.  Loosen/Adjust

12. Entering and Exiting the Contaminated Area

Prior to entering the contaminated area, complete the inspections and checks listed in these user instructions.

1. Ensure Seal Check has been done and motor is responding to your breathing.
2. Check no alarms (battery or filter) are sounding or flashing. If the work environment is noisy you may be unable to hear the alarms. **In this case you must follow the special precautions set out in Section 13.**
3. Don the respirator. **Check other PPE is being worn if required.** Enter the work area
4. Do not remove the respirator until you have left the contaminated area. If you have an acute and spontaneous health episode (ie experiencing dizziness) and believe removing the respirator while you leave the contaminated area may help), remove the respirator with caution.
5. Leave the contaminated area immediately if any of the following conditions occur:
 - Battery or filter alarms trigger.
 - The system is damaged.
 - Airflow in the mask decreases or stops.
 - Breathing becomes difficult
 - You feel dizzy or your vision is impaired.
 - You smell or taste contaminants.
 - Your face, eyes, nose or mouth experience irritation.
 - You suspect the concentration of environmental contaminants has reached levels at which this respirator may no longer provide protection or the oxygen/carbon dioxide levels change so that PAPRs should not be worn.
6. Specific exiting and decontamination procedures may be in place for your workplace as part of the respiratory protection program. Follow those procedures with the instructions below when removing the respirator.

DOFFING THE RESPIRATOR:

Remove the respirator by locating one of the Mask Release Clips (see Section 7) on either side of the respirators. Squeeze the clip and gently pull the respirator from the mask. When the Respirator is removed, the motor/airflow will continue to blow for 10 seconds and then automatically switch to Standby Mode when no breathing is detected. There is no need to turn the respirator off. See Figure 22.


Figure 22 - Mask Release Clips

NOTE: Once In Standby Mode (Green Battery LEDS on), if the respirator is not donned (ie no breath is detected) within three (3) minutes the Respirator will switch into Off Mode (Green LEDS off) and the motor will not be triggered by the wearer's breathing.

DO NOT REMOVE THE RESPIRATOR UNTIL YOU HAVE VACATED THE CONTAMINATED AREA . If you have an acute and spontaneous health episode (ie experiencing dizziness) and believe removing the respirator while you leave the contaminated area may help), remove the respirator with caution.

13. Working in Noisy Environments

CleanSpace respirators have two audible and visual (2) alarms for battery and filters. If the environment is noisy and the user will not hear the alarms, the following precautions must be taken:

- Do not enter the contaminated zone unless all three green battery indicator lights are illuminated
- Do not work in a high noise area for more than four (4) hours. At the end of four (4) hours, leave the contaminated zone and check;
 - The filter blocked light is not illuminated; AND
 - All three green battery indicator lights are still illuminated
- If two or fewer green battery indicator lights are illuminated, recharge CleanSpace Respirator until 3 lights are showing again (and the 3rd one has stopped flashing)
- If the filter blocked light is illuminated, change the filter
- Be alert for breathing difficulty, increased resistance or air flow stopping. If these things happen, exit the contaminated zone immediately.

14. Fitting & Changing the Filter

Before using the respirator you must ensure that you are using the correct filter type for the environment you are working in.

IT IS ESSENTIAL THAT THE CORRECT FILTER TYPE IS SELECTED FOR THE CHOSEN APPLICATION.

- Always change the filter outside of the contaminated area and with the respirator removed
- **CleanSpace filters are NOT reusable and must NOT be cleaned. Cleaning the filter may cause damage.**
- Used filters should be disposed of responsibly and treated as non-recyclable hazardous* waste (dependent on the contaminant being filtered).
- Change the filter if the filter is damaged or there is suspected damage.
- Filters should be changed regularly. The frequency of filter replacement depends on several factors: period of use, concentration of airborne contaminants, exertion levels and workplace protocols where hazards present a risk for cross contamination (ie biohazards or asbestos).

BLOCKED FILTER ALARM

CleanSpace respirators have a Filter Blocked Alarm that triggers when the filter needs to be replaced. The Filter alarm signals by red flashing LED lighting up on the keypad and an audible two (2) beeps, repeated every second. When exiting the contaminated area, the audible Filter Alarm can be muted by pressing the power button once. The audible Filter Alarm will resume if i) after 15 minutes the filter has not been changed or ii) the device is switched on and donned again without the filter being changed.

If the blocked filter alarm is triggered (2 beeps, repeated every second, red LED flashes), leave the contaminated area immediately and replace the filter. Operating the respirator when the blocked filter alarm is sounding can cause the airflow to fall below the manufacturer's minimum designed flow, and may result in mask leak and exposure to contaminants leading to sickness or death.

WHEN TO CHANGE THE FILTER

- Blocked Filter Alarm triggers – audible (2 beeps per second) and visual (Red flashing LEDs).
- The Flow Test indicates that the respirator is not able to produce the Minimum Design Flow.
- The filter is heavily soiled or wet.
- Any sign of damage or suspected damage to the filter.
- The filter reaches its expiry date. The expiry date is marked on the filter label.
- Signs of dust or contaminants on the inside surface of the filter.

A respirator with a clean filter will run for much longer than one with a filter that contains dust contaminants. To maximise your work time between battery charges, change the filter often in high dust environments.


DO NOT USE COMPRESSED AIR OR ANY OTHER MECHANICAL METHOD TO CLEAN THE FILTER! FILTERS ARE VERY EASILY DAMAGED BY THE USE OF COMPRESSED AIR OR BY BRUSHING. Misuse of CleanSpace Respirators by cleaning the filter may result in overexposure to contaminants and lead to sickness or death.

STEPS TO CHANGE THE FILTER ON THE CLEANSPEACE HALO

- To change the filter on the CleanSpace Halo, unlock the filter cover by pulling the pin up. The silver pin is located on the left side of the Respirator (see Figure 23).
- The filter will become loose and fall out. Do not touch the contaminated filter. If necessary, tilt the Respirator down so the filter falls out.
- Visually check the area where the filter sits. Ensure there is no dust or liquids on the inside before replacing the filter. If necessary wipe with a clean cloth or non-alcoholic cleaning wipe.
- To fit a new filter, locate the short edge on the right hand side of the filter. Check the orientation of the "CleanSpace Halo" brand is the right way up. Fit the short edge to the inside (see Figure 23) on the right of the Respirator. Swing the filter shut into the closed position. Hold the filter firmly against the body of the respirator with one hand and press the pin down and back in place. The filter should sit securely in place.

Figure 23 – Replacing the filter in the CleanSpace Halo.

15. Battery Information

CleanSpace HALO is not intrinsically safe. Do not use in flammable or explosive atmospheres. Doing so may result in injury or death.

Always correctly use and maintain the internal lithium ion battery packs. Failure to do so may result in fire or explosion or could adversely affect respirator performance and result in injury, sickness or death.

- Do not charge the internal battery with unapproved chargers or in enclosed cabinets without ventilation or near flammable liquids or gasses, or near sources of high heat
- Do not immerse the device in water
- Do not use, charge or store the device outside the recommended temperature limits

ASSESSING STATE OF BATTERY CHARGE

The level of battery charge is indicated by the three (3) Green LEDs on the keypad. **The Respirator should not be worn when on charge.** The Battery Indicator Lights allow the wearer to estimate the charging level of the battery (see table below).

For the battery to be 100% charged plug in the charger until all three (3) Green LEDs are on and solid (no flashing).

With the unit NOT on charge		
Green LEDs	Charge	Approximate Run Time*
●●●	85 - 100%	4 – 8 hours
●●	15 – 85%	1 – 4 hours
●	5 – 15%	Recharge! (20min to 1 hour)


* *Operating time is affected by mask seal, work rate and other factors. The operating times quoted above are average durations at moderate work rates at sea level. Actual operating times may vary widely from quoted average durations.*

CHARGING THE BATTERY

The battery must only be charged with a CleanSpace Battery Charger.

- To charge the internal battery, turn the respirator upside down.
- Locate the black charging port. See Figure 24 - Charging Port
- Insert the charger connector into the charging port.
- Ensure the charger is plugged into a power source.
- The Battery Indicator Green LEDs on the keypad should be lit up and flashing as the Respirator charges.
- **The Respirator can be charged at temperatures between 0°C to 35°C. The battery will not charge outside this temperature range.**

Figure 24 - Charging Port

Check that at least one of the Battery Indicator Green LEDs on the keypad is on and flashing. If none of the Green LEDs are flashing, the unit is not charging. Check the charger is connected into the Respirator and the power source. If the respirator is still not charging, contact CleanSpace Customer Service sales@cleanspacetechnology.com or visit our website:www.cleanspacehealth.com

- **Charging is complete when the respirator shows three (3) green solid lights.** If the Respirator is unplugged when the third Green LED is flashing, charging is 95% complete.
- When the respirator is fully charged, disconnect the charger cable from the Respirator charging port.
- **IMPORTANT: The Respirator should NOT be worn when on charge.**

Do not charge the battery in electro-magnetic environments as it may damage the electronics.

LOW BATTERY ALARM

CleanSpace respirators have a Low Battery Alarm that is both audible (3 beeps, repeated every second) and visual (Green LEDs) that triggers when the battery has low battery charge. On moderate work rates, this may be 5 – 15 minutes run time (depending on filter loading). Unlike the filter alarm, the Low Battery Alarm cannot be muted. The Alarm will turn off when the respirator is charging or by switching the Respirator off. When the battery charge is extremely low, the motor and airflow will stop.

If the battery alarm sounds (3 beeps, repeated every second), leave the contaminated area immediately and re-charge the battery. Operating the respirator after the low battery alarm has sounded can cause the flow to fall below the manufacturer's minimum designed flow, which may result in mask leak and exposure to contaminants that may lead to sickness or death.

GENERAL BATTERY INFORMATION

CleanSpace respirators have an internal lithium ion (Li-ion) polymer battery. Lithium ion (Li-ion) batteries have the highest energy density of all battery types and are widely used in portable electronic devices. The commonly available Li-ion polymer batteries are packaged in a soft polymer film. The soft packaging allows the cells to expand under extreme heat, making the battery safer in fault conditions. CleanSpace respirators use quality Li-ion polymer batteries and are designed to be viable over a minimum of 500 complete charge and discharge cycles while still maintaining at least 70% of their specified full capacity.

RESPIRATOR STORAGE

To maintain the performance integrity of the internal battery, the respirator must be stored under the following conditions:

Long term storage (> 30 days)	30% - 50% RH, 18°C to 28°C
Short term storage (< 30 days)	30% - 50% RH, -10°C to 35°C

NOTE: Precautions for use:

- Avoid mechanical shocks or impacts from any sharp or hard objects.
- Do not charge or use the respirator if there are any signs of severe mechanical damage
- The respirator will stop functioning if the internal temperature rises above 60°C or falls below -10°C
- Do not dispose of the respirator in a fire.
- Do not use the respirator in heavy rain, or submerged in water or liquids
- Do not disassemble the Respirator case. Disassembly will void the Manufacturer’s Warranty. There are no user serviceable parts inside the Respirator.

IN THE EXTREMELY RARE CIRCUMSTANCE THAT THE BATTERY IS DAMAGED AND ELECTROLYTE COMES IN CONTACT WITH EYES, FLUSH WITH WATER IMMEDIATELY AND SEEK URGENT MEDICAL ATTENTION


16. Cleaning

We recommended that you clean your respirator after every use. The Respirator, mask, neck support and head harness should be cleaned separately. To start, disassemble the mask, neck support and harness from the Respirator before cleaning. **Note: The methods to clean the equipment (outlined below) are not a disinfection or sterilisation protocol. If disinfection or sterilisation is required, contact CleanSpace Customer Service sales@cleanspacetechnology.com or visit our website:www.cleanspacehealth.com and download the recommended protocols for disinfection and sterilization.**

CLEANING THE MASK, HEAD HARNESS AND NECK SUPPORT

Before washing the mask, remove the valve cover over exhalation valve leaf by locating the small inlet located under the lower, inside edge of the valve cover. Use your finger nail to unclip and push up the valve cover. Removing the valve cover allows you to see the valve and check there is no dust or contaminants on the valve. See Figure 25.

Figure 25 - Removing the Mask Exhalation Valve cover

1. Non-Alcohol Disinfectant or Cleaning Wipes

- Non-alcoholic cleaning wipes are impregnated with benzalkonium chloride, a bactericidal solution ideally suited for face pieces that are going to be worn in the near term/same day.
- **Wipes containing alcohol or more stringent cleaning agents should be avoided if the mask is not going to be washed before being worn.**
- Wipe the mask, head harness and neck support inside and out.
- For the mask exhalation valve, with the valve cover removed (using the instructions above), gently lift and wipe the leaf edge. When finished Replace the valve cover before use by snapping the cover back in place.

2. Warm soapy water or an industrial washer

- The mask, head harness and neck support can be washed with a mild detergent in warm water or an industrial washer.
- **Use water and drying temperatures less than 50°C.**
- For hand washing use a soft brush or sponge to remove stubborn dirt or grit.
- For the mask exhalation valve, with the valve cover removed (using the instructions above), gently wash the warm water solution through the valve and use a sponge to gently clean the valve surface
- Rinse the mask, head harness and neck support thoroughly in fresh water to remove cleaning residue.
IMPORTANT: Not rinsing the mask thoroughly, may result in residue from cleaning solution causing irritation to the wearer's skin or may cause the exhalation valve to stick.
- Replace the mask exhalation valve cover before use.

Drying the mask

- After cleaning, allow the mask, head harness and neck support to air dry in a clean environment away from sunlight. Position the mask with the valve facing up to prevent water pooling.
- Do not dry the mask, head harness and neck support by exposing directly to heat i.e. hair dryers or heaters. The mask can also be hand-dried with a clean, lint-free cloth.
- **IMPORTANT: Cleaning can cause the exhalation valve leaf to stick. Before use, check that the exhalation valve leaf operates freely by gently lifting the valve. Replace the mask exhalation valve cover before use.**

WARNING: NEVER DRY THE MASK OR EXHALATION VALVE WITH A CLOTH THAT MAY LEAVE BEHIND LINT. Lint contamination of the exhalation valve may cause it to leak, resulting in overexposure to contaminants and lead to sickness or death.

CLEANING THE RESPIRATOR

- Removing the mask, head harness and neck support. Remove the used filter.
- **IMPORTANT:** Insert the CleanSpace Halo Cleaning & Storage Plug (CS3011) into the filter inlet and two bellow air outlets. Using the Cleaning & Storage Plug (CS3011) prevents dust or liquids from entering the air path of the respirator.
- Cleaning wipes or a cloth soaked with cleaning agent, can be used to wipe down the outside of the respirator. After cleaning, use a clean cloth to remove the cleaning agent.
- Allow the respirator to air dry in a clean environment. The respirator can also be dried with a clean lint-free cloth. Note: **Do not dry the respirator by exposing to heat (i.e. hair dryers or heaters) or compressed air.**
- **When dry, store the Respirator with the Cleaning Storage Plug (CS3011) in place to prevent dust and liquids entering the air path of the respirator.**
- **If the respirator is to be used immediately after cleaning, insert a clean filter and lock in place ready for use.**

WARNING: The respirator contains a battery, sensitive electronics and a motor. NEVER immerse it in water or use anything wetter than a damp cloth to clean it.

WARNING: NEVER CLEAN THE FILTER WITH COMPRESSED AIR OR A BRUSH. FILTERS ARE EASILY DAMAGED BY THE USE OF WATER, COMPRESSED AIR OR BY BRUSHING. Misuse of CleanSpace Respirators by cleaning the filter may result in overexposure to contaminants and lead to sickness or death.

17. Periodic Maintenance & Checking

Recalibration

The respirator can be re-calibrated (see Section 10).

You must recalibrate the internal pressure sensor any time that your CleanSpace Respirator is exposed to SIGNIFICANT CHANGES in temperature of more than 36°F (20°C).). It is best practice to also re-calibrate if the unit has been in storage, particularly if the storage temperature is not known.

18. Approvals

NOTIFIED BODY

Respiratory Protection:

BSI (0086)

Davy Avenue

Knowlhill, Milton Keynes, MK5 8PP

United Kingdom

19. Specifications

Air Flow	CleanSpace Respirators are breath-responsive. Maximum airflow to mask 120 – 220 litre per minute, depending on filter loading, altitude and battery condition.
Operating: Temperatures Humidity Range	Operating Temperature Range: -10°C to 45°C. <i>The respirator will stop if temperature of the unit is greater than 60°C.</i> Humidity Range: -10°C Zero to 90%, non-condensing
Operating: Altitude Range	Approximately sea level to 1770m
Charging: Temperatures	Charging Temperature Range: 0°C to 35°C <i>The battery will not accept charge outside of this range.</i>
Storage Conditions	Short term storage under 30 days: -10°C to 35°C. (30% - 50% RH) Long term storage (> 30 days): 18°C to 28°C (30% - 50% RH) <i>Out of direct sunlight, in a clean, dry environment</i>
Shelf Life	Respirator, mask, head harness and neck support: 3 years (under recommended storage conditions) HEPA Filter: 3 years
Battery Type	Lithium Ion Polymer battery
Battery Charger	Input: 100 – 240VAC, 50 – 60 Hz
Run Time	Up to 8 hours of operating time. <i>Operating time is impacted by work rate, mask leak, filter dust loading, altitude and other factors. Run times may vary widely. IMPORTANT: Short run times (under 4 hours) may be a sign of poor mask fit – seek advice from a trained safety specialist or CleanSpace Support www.cleanspacehealth.com</i>
Recharge Time	2 hours (to 95%)
Low Battery Alarm	Triggers when the Battery Charge is low and run time is around 5 -15 minutes depending on wok load and filter loading. <i>Audible Battery alarm: 3 beeps, repeated every second, 75dB(A) at ear.</i>
Filter Blocked Alarm	Triggers when the filter loading is high and the filter requires replacement. <i>Audible alarm: 2 beeps repeated every second, 75dB(A) at ear. Plus red filter LED flashes.</i>
Weight	Complete system: 533g (Respirator (CS3001), Filter, Mask, Neck Support, head harness)
Intrinsic Safety	CleanSpace HALO is not an intrinsically safe system

20. Product and Accessory information

Respirators

Product Code	Description
CS3001	Respirator includes case, internal battery, motor and electronics.

Filters

Product Code	Description
CS3002	Filter HEPA/P3/TM3 Particulate Glass

Masks

Product Code	Description
CS3003	Half Mask – Small
CS3004	Half Mask – Medium
CS3005	Half Mask – Large
CS3006	Full Face Mask - Small
CS3007	Full Face Mask – Medium/Large

Accessories – Head Harness, Neck Supports, Charger and Cleaning Plug

Product Code	Description
CS3008	Head Harness
CS3009	Neck Support - Small
CS3010	Neck Support - Medium
PAF-1101	Battery Charger Universal – Multi
PAF-1009	Seal Check Cap
CS3011	Cleaning and decontamination Plug
CS3013	CleanSpace Respirator ID panels (different colours)
CS3015	CleanSpace Mask ID badges (different colours)
CS3014	Charging & Storage Station (Charges and holds 8 units)

21. Product Warranty

This product has been manufactured using quality parts and processes. CleanSpace Technology Pty Ltd warrants that the product is free from defective workmanship and parts for a period of three (3) years from the date of original purchase provided the product has been used, cleaned and maintained in accordance with these instructions and CleanSpace's recommendations. This warranty does not include consumable parts, such as filters and face masks, which must be replaced regularly by the user. Consumable parts are warranted up to the point of use, provided they have been stored correctly and are within their expiry date.

This warranty does not cover:

- Where the product has been used for industrial purposes outside the recommendations of CleanSpace Technology Pty Ltd;
- Where damage has been caused by misuse, neglect, accident, or excessive wear and tear.

Any claim under this warranty must be made within three (3) years of the date of purchase of the product. All warranty claims must be made by returning the defective product to your supplier together with the proof of purchase. The purchaser is responsible for all freight. In the event that any part of the product is found by CleanSpace to be defective, CleanSpace will either repair or at its discretion replace the faulty part.

This warranty is given by:

CleanSpace Technology Pty Ltd (ABN 24 146 453 554),
A Ground Floor, 16-18 Carlotta Street Artarmon NSW 2064 Australia;
T. +61 2 8436 4000 E. support@cleanspacetechnology.com

This warranty is provided in addition to other rights and remedies you have under law. You are entitled to replacement or refund for a major failure. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Disclaimer:

Whilst CleanSpace Technology Pty Ltd has made every effort to ensure that the details and information given in both our printed and online publications are accurate at the time of issue, full technical specifications are not necessarily included. Furthermore CleanSpace Technology Pty Ltd has a policy of continuous improvement and the right is reserved to alter details and information as the need arises. Accordingly the Customer should check any details and information they wish to rely on with CleanSpace Technology Pty Ltd at the time of purchase. CleanSpace Technology Pty Ltd cannot accept liability in respect of any errors or omissions herein contained or for any loss or damage malfunction or consequential loss arising from reliance upon our publication.

The Customer will be responsible for any risk to health or safety from goods in the Customer's possession and/or control. The Customer's attention is drawn to the fact that statutory regulations and recognised codes of practice exist covering the use and handling of some goods (including safety products). The Customer must ensure that persons who use the goods receive adequate training and safety literature.

©2018 CleanSpace Technology Pty Ltd. All rights reserved.

No page or part of this manual may be reproduced in any form without written permission of the copyright owner shown above. CleanSpace reserves the right to correct typographical errors. All information is correct at the time of going to print.

CleanSpace Technology Pty Ltd

ABN 24 146 453 554

Ground Floor, 16-18 Carlotta Street, Artarmon NSW 2064 Australia

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